


The **IntraStent[®] LD** Family
Large Lumen Stents



The **Power of Flexible Strength**

Diverse Selection


ev3 Inc. offers a wide variety of Large Diameter Stents to serve the varied needs of our customers and their patients.

Innovative Design

With the LD series, ev3 engineers have created a product group that combines the features of flexibility, strength and coverage in ways not found in other large lumen stents.

Passion for Results

The people of ev3 are committed to the notion that growth and innovation in response to our customers will lead to consistently superior performance, clinical satisfaction and patient outcomes.





IntraStent® DoubleStrut™ LD Stent

The IntraStent DoubleStrut LD Stent is our original large diameter stent design. Upon its introduction, this device set industry standards for flexibility and profile, leading to superior deliverability and changed the way clinicians viewed large diameter stents.

- Unique DoubleStrut design allows unbeatable flexibility
- Low profile, 8F Introducer sheath
- No shortening after expansion up to 12 mm
- Rounded edges minimize potential for lumen trauma

When flexibility is your greatest clinical concern, try the IntraStent® DoubleStrut™ LD Stent






IntraStent® Mega™ LD Stent

The second offering in the IntraStent LD Family is the Mega LD Stent. This stent combines the best features of the DoubleStrut LD Stent and the Max LD Stent to create a device that is strong, yet pliant.

- 5 cell design for increased radial strength over the Doublestrut LD Stent
- 9F sheath compatibility
- No shortening after expansion up to 12 mm
- Rounded edges minimize potential for lumen trauma

When the clinical situation demands strength AND flexibility, look to the IntraStent Mega LD Stent



IntraStent® Max™ LD Stent

The third offering in the IntraStent LD Family is the Max LD Stent. This stent is engineered for strength and coverage, yet is still very deliverable.

- 6 cell design for maximum strength and coverage
- 11 F sheath compatibility
- No shortening after expansion to 12 mm
- Rounded edges minimize potential for lumen trauma

When compressive forces are putting the clinical situation at risk, reach for the IntraStent Max LD Stent

	IntraStent® DoubleStrut™ LD	IntraStent® MEGA™ LD	IntraStent® MAX™ LD
Strength	★	★★	★★★
Flex	★★★	★★	★★
Shortening	None 9-12mm	None 9-12mm	None 12mm
Profile	★★★ 8Fr	★★ 9 Fr	★ 11 Fr

Ordering Information

Product Number	Unexpanded		Expanded	
	Stent Dimensions Diameter (mm)	Stent Dimensions Length (mm)	Stent Dimensions Diameter (mm)	Stent Dimensions Length (mm)
IntraStent® DoubleStrut™ LD Stent				
S15-16	3.8	16.0	9,10,11,12	16.0
S15-26	3.8	26.0	9,10,11,12	26.0
S15-36	3.8	36.0	9,10,11,12	36.0
S15-56	3.8	56.0	9,10,11,12	56.0
S15-76	3.8	76.0	9,10,11,12	76.0
IntraStent® Mega™ LD Stent				
S17-16	3.8	16.0	9,10,11,12	16.0
S17-26	3.8	26.0	9,10,11,12	26.0
S17-36	3.8	36.0	9,10,11,12	36.0
IntraStent® Max™ LD Stent				
S18-16	4.5	16.0	12	16.0
S18-26	4.5	26.0	12	26.0
S18-36	4.5	36.0	12	36.0

Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

These stents are CE marked for use in iliac and subclavian arteries. IntraStent, DoubleStrut, Mega and Max are trademarks of **ev3** Inc.



AN ENDOVASCULAR COMPANY

4600 Nathan Lane North
Plymouth, MN 55442-2920
USA
PH 763.398.7000
FX 763.398.7001
Toll Free: 800.716.6700

ev3 Europe SAS
Paris Nord II
9 rue des 3 Soeurs
BP 60002 Villepinte
95970 ROISSY CDG CEDEX
France
PH +33 (0) 149 89 00 22
FX +33 (0) 149 89 00 33

France
PH +33 (0)149 89 00 22
FX +33 (0)149 89 00 33
Italy
PH +39 0331 781 212
FX +39 0331 780 499

UK
PH +44 1279 659 900
FX +44 1279 654 900

Germany,Austria,Switzerland
PH +49 228 52883 0
FX +49 228 52883 60

Benelux
PH +31 433 659 223
FX +31 433 646 395

Nordic
PH +46 878 325 60
FX +46 878 325 69

Spain
PH +34 91 656 7154
FX +34 91 656 7214

Portugal
PH +35 1 214 170 081
FX +35 1 214 170 080