


Visi-Pro™
**Balloon Expandable
Peripheral Stent System**




Visible *precision.*

Visible Precision.

Visi-Pro™ Balloon Expandable Peripheral Stent System


Premium Visibility

Exclusive 0.035" BE Stent with radiopaque stent markers ensuring visualization during and after stent placement





Precise Placement

- Precise balloon to stent placement minimizing dilatation of healthy tissue
- Accurate catheter marker alignment to balloon taper
- Progressive stent design delivers minimal shortening



Solid Radial Strength and Flexibility

- Patented stent design for strength and flexibility¹



Superior Trackability and Retraction

- 5-Fold Balloon wrap improves retraction
- 28% slimmer catheter tip profile for tracking tortuous anatomy**
- All 8 mm stents and smaller are 6F compatible

¹ Comparison testing performed by ev3 Inc. Data on file. Bench results not intended to indicate clinical performance. **Comparable to Primus™ GPS™

Visi-Pro™ Balloon Expandable Peripheral Stent System

Ordering Information

Product Number Catheter Length 80 cm	Product Number Catheter Length 135 cm	Stent Diameter (mm)	Stent Length (mm)	Balloon Length (mm)	Sheath Size (F)
PXP35-05-12-080	-	5	12	15	6*
PXP35-05-17-080	PXP35-05-17-135	5	17	20	6
PXP35-05-27-080	PXP35-05-27-135	5	27	30	6
PXP35-05-37-080	PXP35-05-37-135	5	37	40	6
PXP35-05-57-080	PXP35-05-57-135	5	57	60	6
PXP35-06-12-080	-	6	12	15	6
PXP35-06-17-080	PXP35-06-17-135	6	17	20	6
PXP35-06-27-080	PXP35-06-27-135	6	27	30	6
PXP35-06-37-080	PXP35-06-37-135	6	37	40	6
PXP35-06-57-080	PXP35-06-57-135	6	57	60	6
PXP35-07-12-080	-	7	12	15	6
PXP35-07-17-080	PXP35-07-17-135	7	17	20	6
PXP35-07-27-080	PXP35-07-27-135	7	27	30	6
PXP35-07-37-080	PXP35-07-37-135	7	37	40	6
PXP35-07-57-080	PXP35-07-57-135	7	57	60	6
PXP35-08-17-080	PXP35-08-17-135	8	17	20	6
PXP35-08-27-080	PXP35-08-27-135	8	27	30	6
PXP35-08-37-080	PXP35-08-37-135	8	37	40	6
PXP35-08-57-080	PXP35-08-57-135	8	57	60	6
PXP35-09-17-080	PXP35-09-17-135	9	17	20	7
PXP35-09-27-080	PXP35-09-27-135	9	27	30	7
PXP35-09-37-080	PXP35-09-37-135	9	37	40	7
PXP35-09-57-080	PXP35-09-57-135	9	57	60	7
PXP35-10-17-080	PXP35-10-17-135	10	17	20	7
PXP35-10-27-080	PXP35-10-27-135	10	27	30	7
PXP35-10-37-080	PXP35-10-37-135	10	37	40	7
PXP35-10-57-080	PXP35-10-57-135	10	57	60	7

*6F=0.085" I.D.

Visi-Pro™ Compliance Chart

Visi-Pro™ Diameter (mm)	Inflation Pressure (atm)				
	8	9	10	11	12
5.0	5.00	5.09	5.16	5.22	5.28
6.0	6.00	6.11	6.22	6.31	6.39
7.0			7.00	7.09	7.17
8.0			8.00	8.15	8.26
9.0			9.00	9.15	9.28
10.0			10.00	10.11	10.21

● Diameter at Nominal Pressure ● Diameter at Rated Burst Pressure

Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

The Visi-Pro™ is CE marked for use in the renal, common and external iliac, and subclavian arteries, as well as malignant biliary use.

Visi-Pro™ is a trademark of ev3 Inc.

NEUROVASCULAR | PERIPHERAL VASCULAR

Access - Balloons - Carotid - Embolic Coils - Embolic Protection - Liquid Embolics - Plaque Excision - Procedural Support - Retrieval Devices - Stents

ev3 Europe

International Headquarters
106-108 rue La Boétie
75008 Paris
France
PH +33 156 88 59 10
FX +33 156 88 59 11

www.ev3.net

ev3 Corporate

World Headquarters
Peripheral Vascular
3033 Campus Drive
Plymouth, MN 55441
USA
PH +1 763 398 7000
FX +1 763 398 7001

ev3 Neurovascular

9775 Toledo Way
Irvine, CA 92618 USA
PH +1 949 837 3700
FX +1 949 837 2044

ev3 International

Distribution Centre
Europalaan 25
6199 AB Maastricht-Airport
The Netherlands
PH +31 (0) 433 659 220
FX +31 (0) 43 364 6395

ev3 SAS France

PH +33 (0) 156 88 31 10
FX +33 (0) 156 88 31 11

ev3 B.V. Benelux

PH +31 (0) 433 659 223
FX +31 (0) 433 650 283

ev3 Technologies Iberica, S.L. Spain

PH +34 91 656 7154
FX +34 91 656 7214

ev3 S.r.l Italy

PH +39 0267 977 61
FX +39 0266 711 637

ev3 Nordic AB

PH +46 859 000 950
FX +46 859 000 959

ev3 Sp z o.o. Poland

PH +48 32 747 01 44
FX +48 32 747 01 45

ev3 GmbH Germany, Austria

PH +49 228 528 830
FX +49 228 528 8360

ev3 Ltd. United Kingdom

PH +44 1279 659 900
FX +44 1279 654 900

EC REP
MediMarkt Europe Sarl
BP 2332
F-38033 Grenoble Cedex 2
FRANCE

